

# Using $\mu$ Manager supported devices in Zen

Jérôme Mutterer<sup>1\*</sup>, Julien Kissenberger<sup>2</sup> and Fabrice Schmitt<sup>2</sup>

<sup>1</sup>C.N.R.S., Strasbourg, France.

<sup>2</sup>Carl Zeiss S.A.S., Marly-le-Roi, France.

\*Corresponding author(s). E-mail(s): [jerome.mutterer@cnrs.fr](mailto:jerome.mutterer@cnrs.fr);  
Contributing authors: [julien.kissenberger@zeiss.com](mailto:julien.kissenberger@zeiss.com);  
[fabrice.schmitt@zeiss.com](mailto:fabrice.schmitt@zeiss.com);

## Abstract

Microscopy setups can include common or custom devices for which no driver is available in the Zen software. In this white paper, we propose a method and show an example of creating a basic graphical user interface for driving those peripherals from Zen, provided a matching  $\mu$ Manager device adapter is available.

**Keywords:** drivers, adapters, devices,  $\mu$ Manager, pycromanager, Zen

## 1 Introduction

Advanced microscopy experiments conducted in core imaging facilities typically require systems of two very distinct types. First, we find complete commercially available systems. Those usually consist of official manufacturer hardware parts like the microscope stand, beam path mirrors, objectives or reflector turrets. Motorized versions of those core elements are of course well integrated into the company's provided microscope control software. For advanced experiments, these core elements can be complemented by so-called third party provided pieces of hardware, as is often the case for motorized stage, specific illuminators, additional filter wheels, or environmental control devices. The ability of the main control software to communicate with those third party devices will make user training and workflows execution easier, as a single software can be used throughout the experiment.

Home-made, or custom microscopy setups make up the second type of core facilities imaging systems. Those can be anything ranging from 3D printed structural parts to a mix of up cycled outdated stands, standalone filter wheels and Arduino driven illuminators or motion controllers.  $\mu$ Manager, the Open Source microscope control and automation software [1, 2] is often used to try and integrate those diverse hardware parts, given its extensive list of so called device adapters supporting direct communication with the devices themselves or talking to the system level device drivers.

While the current version of  $\mu$ Manager is impressively capable for advanced acquisition protocols, commercially available control software like Zen [3] still benefit from smoother user experience, with greater development efforts targeted to user experience.

Here we present a simple step-by-step method to create a Zen graphical user interface that drives devices configured to work in  $\mu$ Manager, enabling creating imaging setups comprising of both vendor supported and  $\mu$ Manager supported hardware, using the `pycromanager`[4] Python library.

## 2 Overview

1. Prerequisite: Zen 3.4 (blue edition), with Macro Environment Module
2. Install  $\mu$ Manager and configure devices
3. Install Python 3.6
4. Install and test `pycromanager`
5. Create Python scripts to query and alter device state
6. Create Zen Macro that invokes scripts from step 5

## 3 Walk-through

### 3.1 Prerequisite: Zen 3.4 (blue edition), with Macro Environment Module

This tutorial is intended for a working imaging setup running Zen 3.4. It uses Zen OAD Macros, so make sure the Macro Environment Module is activated. To do this, navigate to Zen>Tools>Modules Manager...




Fig. 1 Zen optional modules manager

### 3.2 Install $\mu$ Manager and configure devices

In this section, we install  $\mu$ Manager and configure a sample device according to instructions available from the  $\mu$ Manager wiki at <https://micro-manager.org/>

Retrieve the latest  $\mu$ Manager software from the [download page](#). For this tutorial, we chose the Windows 64-bits version [version 2.0.0](#).

Click through the installer steps and launch  $\mu$ Manager.

In the "Micro-Manager Startup Configuration" dialog, select "Default User" profile and "(none)" as "Hardware Configuration File" as we're going to add just the devices we need.





Fig. 2  $\mu$ M Startup Configuration Dialog

$\mu$ Manager starts with no devices configured.



**Fig. 3** Main  $\mu$ M user interface

The example device we will use is a custom made device emulating a PreciExcite LED illuminator. It is built around an Arduino micro-controller and a NeoPixels RGB LED ring. It is connected to the host computer using a serial over USB port, and the micro-controller firmware implements a subset of the official PreciExcite protocol. This allows  $\mu$ Manager hardware configuration tool to recognize the device as a 3 channels illumination device with shutter. Further examples of emulated custom devices and their Arduino firmwares can be found at [https://github.com/mutterer/MM\\_CustomArduino](https://github.com/mutterer/MM_CustomArduino)




**Fig. 4** A custom illumination device built around an Arduino micro-controller

Now start the configuration assistant by navigating to Micro-Manager>Devices>Hardware Configuration Wizard...

1. Chose "Create new configuration" and press "Next"
2. From the devices list, open the "PreciExcite" folder, select the "PreciExcite: PreciExcite LED illuminator" device and add it using the "Add..." button



The extensive documentation at [pymanager](http://pymanager.org) website provides a test program to make sure it can communicate with  $\mu$ Manager. Open a Python command line window, breath, and type the following code:



```
Python 3.6 (64-bit)
Python 3.6.0 (v3.6.0:41df79263a11, Dec 23 2016, 08:06:12) [MSC v.1900 64 bit (AMD64)] on win32
Type "help", "copyright", "credits" or "license()" for more information.
>>> from pymanager import Bridge
>>> with Bridge() as bridge:
...     print(bridge.get_core())
...     print(bridge.get_core().get_shutter_device())
...
<pymanager.core.mmcorej_CMMCore object at 0x000001F353060DA0>
PrecisExcite
>>>
```

**Fig. 6** Testing pymanager successful installation


We were able to programmatically get a reference to  $\mu$ Manager core and from this the name of the current shutter device we installed earlier.

## 3.4 Python scripts

It is unfortunate that we cannot directly use pymanager from Zen macros, so we first need to create two intermediate scripts.

### 3.4.1 Device query script

The first one queries the current illuminator states, the 3 channels intensities and the shutter state and replies with a simple text string containing those values. The following code is basic and can easlily be improved.



```
test02_readSettings.py - C:\Users\jmutterer\Desktop\test02_readSettings.py (3.6.0)
File Edit Format Run Options Window Help
import sys
from pymanager import Bridge
with Bridge() as bridge:
    # get a reference to MM core
    core = bridge.get_core()
    # read current device properties
    b=core.get_property("PrecisExcite","IntensityA")
    g=core.get_property("PrecisExcite","IntensityB")
    r=core.get_property("PrecisExcite","IntensityC")
    s=core.get_property("PrecisExcite","State")
    # output properties as a semicolon separated string
    print (b,g,r,s,sep=";")
Ln: 13 Col: 0
```

**Fig. 7** A python script that queries device state and reports it


You can test this script by opening a Windows command line and typing:

```
python Desktop/test02_readSettings.py
10;20;8;1
```

The replied string means intensities for channels A,B,C are 10,20,8, and the final 1 means the shutter is open (LEDs are shining).

### 3.4.2 Device driving script

The second script is in charge of setting device properties. It must be called with arguments and uses pycromanager to send the matching commands to  $\mu$ Manager.



```

test02_applySettings.py - C:\Users\jmutterer\Desktop\test02_applySettings.py (3.6.0)
File Edit Format Run Options Window Help
import sys
from pycromanager import Bridge
with Bridge() as bridge:
    # get a reference to MM core
    core = bridge.get_core()
    # set device properties as per recieved script arguments
    core.set_property("PrecisExcite", "IntensityA", sys.argv[1])
    core.set_property("PrecisExcite", "IntensityB", sys.argv[2])
    core.set_property("PrecisExcite", "IntensityC", sys.argv[3])
    core.set_property("PrecisExcite", "State", sys.argv[4])
    # optional: refresh MM UI
    studio = bridge.get_studio()
    app = studio.get_application()
    app.refresh_gui()
Ln: 11 Col: 29

```

Fig. 8 A python script that sets device states from arguments it receives

You can test this script by opening a Windows command line and typing:

```
python Desktop/test02_applySettings.py 50 0 50 1
```

There's no reply for this script, but the device should display a bright purple color with 50% blue, 0% green, 50% red, and an open shutter state (1). Congratulations, we have all pieces to put everything together in a Zen macro!

## 3.5 Zen macro

### 3.5.1 A basic example


The last part to design is a Zen macro that will:

1. Execute the first script from 3.4.1 to query the device state
2. Build a user interface with sliders for channels intensities and a shutter checkbox
3. Execute the second script from 3.4.2 to alter the device state

Open Zen. Navigate to Zen>Macro>Macro Editor... and enter the following code:

We kept code to a minimum on purpose, but it will be easy to extend this model.

In a basic setting, you could easily create a Zen Toolbar button item to run this code by using the Macro Editor's "Toolbar Configuration" widget and associate a new toolbar button with the current macro.



**Fig. 9** The Zen macro wrapping calls to pycromanager scripts and displaying a user interface

### 3.5.2 Automation use case

For more complex experiment workflows, a first approach could be to use the Zen built-in "Automation" module. This module allows to execute given macros before and after experiments. A use case would be to pause fluid flow in a perfusion imaging chamber, driving an  $\mu$ Manager configured peristaltic pump.




**Fig. 10** The "Automation" dialog specifying pre- and post-experiment actions (macros)

### 3.5.3 Experiment feedback use case

Whenever actions should be taken during experiment, driving external devices could be integrated in the "Experiment Feedback" workflow design. Experiment feedback provides acquisition hack points at experiment startup and completion, but also allows Zen Macro code to be be upon events occurring at the "observable" level. Observables are system or environment properties whose values are monitored and can be used to trigger code execution. The following example will use those acquisition workflow entry points to adjust illumination level for each time point in an sample experiment and finish experiment by resetting the device after a given value was reached.

First in your "Acquisition" tab, enable "Time Series" and "Experiment Feedback" check boxes. In the matching dialogs, Set time series duration to "As long as possible". Set "Script runtime condition" to "Free Run".



**Fig. 11** Setting up a time lapse experiment with feedback

Now, open the "Feedback Script Editor" and enter the following pre-loop, loop, and post-loop scripts.

The PreLoop script initializes an intensity variable and paths to the Python 3.6 interpreter and pycromanager script responsible to drive our device through  $\mu$ Manager.

The Loop script monitors the "CurrentTimePointIndex" observable and runs the pycromanager script whenever the observable changes. It also increments the intensity, and stops the acquisition when a given threshold is crossed.

Finally, the PostLoop script turns off our device after the acquisition.

```

Script Editor for Experiment Feedback
PreLoop Script - Single Execution on Experiment Start
1 intensity = 0
2 pycro_script = "C:/Users/jmutterer/Desktop/test02_applySettings.py"
3 python3 = "C:/Users/jmutterer/AppData/Local/Programs/Python/Python36/python.exe"
4 ZenService.Xtra.System.WriteDebugOutput("Experiment started")
5

Loop Script - Repetitive execution during experiment runtime whenever an used observable has changed
1 if ZenService.Experiment.HasChanged("CurrentTimePointIndex"):
2     ZenService.Xtra.System.WriteDebugOutput("setting intensity to: "+ str(intensity))
3     ZenService.Xtra.System.ExecuteExternalProgramBlocking(python3, pycro_script+" "+str(intensity)+" 0 0 1",3000)
4     intensity = intensity + 10
5     if intensity > 100:
6         ZenService.Actions.StopExperiment()

PostLoop Script - Single Execution on Experiment Stop
1 ZenService.Xtra.System.ExecuteExternalProgram(python3, pycro_script+" 0 0 1")
2 ZenService.Xtra.System.WriteDebugOutput("Experiment finished")

```

**Fig. 12** Experiment Feedback scripts

## 4 Conclusion

Not all commercially available and basically none of custom made devices benefit from integration with Zen software. In this white paper, we present a way to drive devices from Zen, provided they can be configured to work with  $\mu$ Manager in the first place. The extensive list of supported devices in  $\mu$ Manager and the ability to easily create  $\mu$ Manager device adapters, or emulate existing similar devices opens opportunities for using a wide range of such devices in Zen. The proposed method uses OAD Zen macros and interface python scripts using the pycromanager library to talk to  $\mu$ Manager, and illustrate how Open Source and scriptable vendor software can cooperate to achieve better device support and improved user experience.

## References

- [1] Edelstein, A., Amodaj, N., Hoover, K., Vale, R., Stuurman, N.: Computer control of microscopes using  $\mu$ manager. *Current Protocols in Molecular Biology* **92**(1), 14–201142017 (2010). <https://doi.org/10.1002/0471142727.mb1420s92>
- [2] Edelstein, A.D., Tsuchida, M.A., Amodaj, N., Pinkard, H., Vale, R.D., Stuurman, N.: Advanced methods of microscope control using manager software. *Journal of Biological Methods* **1**(2), 10 (2014). <https://doi.org/10.14440/jbm.2014.36>
- [3] Zen 3.4. Zeiss Efficient Navigation, Zeiss Microscopy (2021). <https://www.zeiss.com/microscopy/int/products/microscope-software/zen.html>
- [4] Pinkard, H., Stuurman, N., Ivanov, I., Anthony, N., Ouyang, W., Li, B., Yang, B., Tsuchida, M., Chhun, B., Zhang, G., Mei, R., Anderson, M., Shepherd, D., Hunt-Isaak, I., Dunn, R., Jahr, W., Kato, S., Royer, L., Thiagarajah, J., Eliceiri, K., Lundberg, E., Mehta, S., Waller, L.: Pycromanager: open-source software for customized and reproducible microscope control. *Nature Methods* **18**(3), 226–228 (2021). <https://doi.org/10.1038/s41592-021-01087-6>